



Human FLT1 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6096
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Ile328
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	37.77 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P17948
ALIAS/SYNONYMS	FLT1, Vascular endothelial Growth Factor proteins receptor 1, VEGFR1, Vascular permeability factor receptor, FRT, Tyrosine-protein kinase receptor FLT, VEGFR-1, FLT, FLT-1, Tyrosine-protein kinase FRT, Fms-like tyrosine kinase 1, Drug Target

