



Human EPHA3 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6099
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Gln541
TAG	C-Terminal Fc Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	87.36 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P29320
ALIAS/SYNONYMS	HEK, Eph-like tyrosine kinase 1, EK4, EPHA3, Tyrosine-protein kinase receptor ETK1, Tyrosine-protein kinase TYRO4, Ephrin type-A receptor 3, hEK4, ETK, TYRO4, ETK1, Human embryo kinase, EPH-like kinase 4, Drug Target

