



Human AXL-UFO recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6100
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Ala26-Pro449
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	49.62 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P30530
ALIAS/SYNONYMS	UFO, AXL oncogene, Tyrosine-protein kinase receptor UFO, AXL, Drug Target, SARS-like Coronaviruses

