



Human THPO recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6102
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Ser22-Gly353
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	38.74 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P40225
ALIAS/SYNONYMS	Myeloproliferative leukemia virus oncogene ligand, C-mpl ligand, Thrombopoietin, THPO, Megakaryocyte colony-stimulating factor, Megakaryocyte growth and development factor, ML, MGDF

