



Human CDw210a-IL10RA recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6105
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Asn235
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	27.51 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q13651
ALIAS/SYNONYMS	IL-10 receptor subunit alpha, IL-10R1, Interleukin-10 receptor subunit alpha, Interleukin-10 receptor subunit 1, CD210, IL10RA, CDw210a, IL-10R subunit alpha, IL-10R subunit 1, IL10R, IL-10RA

