



Human CD85k-LILRB4/ILT3 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6107
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Glu259
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	29.45 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q8NHJ6
ALIAS/SYNONYMS	CD85 antigen-like family member K, LIR-5, Leukocyte immunoglobulin-like receptor 5, ILT3, ILT-3, Immunoglobulin-like transcript 3, LILRB4, CD85k, Leukocyte immunoglobulin-like receptor subfamily B member 4, Monocyte inhibitory receptor HM18, LIR5

