



Human NECTIN4-PVRL4 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6108
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Ser349
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	38.41 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q96NY8
ALIAS/SYNONYMS	LNIR, NECTIN4, PRR4, Poliovirus receptor-related protein 4, Nectin cell adhesion molecule 4, Ig superfamily receptor LNIR, Nectin-4, PVRL4, Drug Target

