



Human CD360-IL21R recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6110
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Glu232
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	27.56 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q9HBE5
ALIAS/SYNONYMS	IL21R, NILR, Interleukin-21 receptor, IL-21R, IL-21 receptor, CD360, Novel interleukin receptor

