



Human IL1RAP recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6111
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Glu359
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	42.38 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q9NPH3
ALIAS/SYNONYMS	IL-1R3, C3orf13, Interleukin-1 receptor 3, Interleukin-1 receptor accessory protein, IL1RAP, IL-1R-3, IL-1 receptor accessory protein, IL1R3, IL-1RAcP, Drug Target

