



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Human IL35 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6112
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	(Arg23-Ser219) & (Arg21-Lys229)
TAG	C-Terminal Fc Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	74.41 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P29459 & Q14213
ALIAS/SYNONYMS	Interleukin-35, IL35, IL-35



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.