



Human SCARB1-SR-BI/CLA-1/SRB1 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6114
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Pro33-Tyr443
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	49.63 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q8WTV0
ALIAS/SYNONYMS	CD36, CLA1, CD36 antigen-like 1, CD36L1, SRB1, Collagen type I receptor, thrombospondin receptor-like 1, SCARB1, Scavenger receptor class B member 1, CLA-1, SR-BI, CD36 and LIMPII analogous 1, SARS-like Coronaviruses

