



Human PADI2 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6117
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Pro665
TAG	N-Terminal Strep Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	78.79 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q9Y2J8
ALIAS/SYNONYMS	KIAA0994, Peptidylarginine deiminase II, PADI2, PAD-H19, PDI2, Protein-arginine deiminase type II, Protein-arginine deiminase type-2, PAD2

