



Human PIP4K2B-PIP4KII-beta recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6121
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Ser2-Thr416
TAG	N-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	48.42KD
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P78356
ALIAS/SYNONYMS	1-phosphatidylinositol 5-phosphate 4-kinase 2-beta, Phosphatidylinositol 5-phosphate 4-kinase type-2 beta, PtdIns(5)P-4-kinase isoform 2-beta, Phosphatidylinositol 5-phosphate 4-kinase type II beta, PIP4KII-beta, Diphosphoinositide kinase 2-beta, PI(5)P 4

