



Human DSG1 recombinant protein (in Insect)

OVERVIEW

PRODUCT REFERENCE	PX-P6123
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Glu50-Pro548
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Insect
MOLECULAR WEIGHT	58.34 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q02413
ALIAS/SYNONYMS	DGI, Cadherin family member 4, Desmoglein-1, DG1, Pemphigus foliaceus antigen, CDHF4, DSG1, Desmosomal glycoprotein 1

