



Monkey CD49b & CD61- ITGA2 & ITGB3 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6125
ORIGIN SPECIES	Monkey
FRAGMENT	Met1-Pro908(CD49b) & Met1-Pro717(CD61)
TAG	C-Terminal His & C-Strep Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	99.55 kDa,
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	EHH58199.1 & XP_014974842.2
ALIAS/SYNONYMS	Collagen receptor, ITGA2, VLA-2 subunit alpha, CD49B, CD49 antigen-like family member B, Platelet membrane glycoprotein Ia, GPIa, Integrin alpha-2, CD49b, GPIIIa, Integrin beta-3, CD61, ITGB3, Platelet membrane glycoprotein IIIa, GP3A

