



Human LGR4 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6128
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Thr544
TAG	C-Terminal Fc Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	85.79 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q9BXB1
ALIAS/SYNONYMS	GPR48, Leucine-rich repeat-containing G-protein coupled receptor 4, LGR4, G-protein coupled receptor 48, Drug Target

