



Human C1S recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6130
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Asp688
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	77.66 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P09871
ALIAS/SYNONYMS	C1 esterase, Complement C1s subcomponent, C1S, Complement component 1 subcomponent s, Drug Target

