



Human LRRC15 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6131
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Gly538
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	60.99 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q8TF66
ALIAS/SYNONYMS	Leucine-rich repeat protein induced by beta-amyloid homolog, LRRC15, Leucine-rich repeat-containing protein 15, LIB, hLib, Drug Target

