



African swine fever virus ASFV E184R recombinant protein

OVERVIEW

| | |
|---------------------------|---------------------------|
| PRODUCT REFERENCE | PX-P6136 |
| ORIGIN SPECIES | African swine fever virus |
| FRAGMENT | Met1-Glu184 |
| TAG | C-Terminal His Tag |
| DELIVERY CONDITION | Dry Ice |

PRODUCT INFORMATION

| | |
|--------------------------------|--|
| EXPRESSION SYSTEM | Eukaryotic expression |
| HOST | Insect |
| MOLECULAR WEIGHT | 22.72kDa |
| FORM | |
| BUFFER | 0.01M PBS, pH 7.4. |
| STABILITY & STORAGE | 4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection) |
| UNIPROT | YP_009927248.1 |
| ALIAS/SYNONYMS | hypothetical protein IM014_gp157, E184L CDS, African Swine Fever Virus (ASFV) |

