



African swine fever virus ASFV F317L recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6137
ORIGIN SPECIES	African swine fever virus
FRAGMENT	Met1-Thr317
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Insect
MOLECULAR WEIGHT	37.55kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	QPB67569.1
ALIAS/SYNONYMS	F317L, African Swine Fever Virus (ASFV)

