



## Human LOXL2 recombinant protein

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-P6142
<b>ORIGIN SPECIES</b>	Homo sapiens (Human)
<b>FRAGMENT</b>	Met1-Gln774
<b>TAG</b>	C-Terminal His Tag
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Eukaryotic expression
<b>HOST</b>	Mammalian cells
<b>MOLECULAR WEIGHT</b>	87.69 kDa
<b>FORM</b>	
<b>BUFFER</b>	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
<b>STABILITY &amp; STORAGE</b>	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
<b>UNIPROT</b>	Q9Y4K0
<b>ALIAS/SYNONYMS</b>	LOXL2, Lysyl oxidase-related protein WS9-14, Lysyl oxidase homolog 2, Lysyl oxidase-like protein 2, Lysyl oxidase-related protein 2

