



Human MGAM-MGA recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6144
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Asp960-Trp1853
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	64.21 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	O43451-1
ALIAS/SYNONYMS	Maltase-glucoamylase, intestinal, Glucan 1,4-alpha-glucosidase, MGA, Alpha-glucosidase, MGAML, MGAM

