



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Human FGF9 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6145
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Leu4-Ser208
TAG	N-Terminal Fc Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	50.97 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P31371
ALIAS/SYNONYMS	FGF9, Glia-activating factor, GAF, Heparin-binding Growth Factor proteins 9, FGF-9, Fibroblast Growth Factor proteins 9, HBGF-9



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.