



Human BMP4 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6146
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Arg292-Arg408
TAG	N-Terminal Fc Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	40.96 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P12644
ALIAS/SYNONYMS	BMP2B, BMP-2B, DVR4, Bone morphogenetic protein 4, BMP-4, BMP4, Bone morphogenetic protein 2B

