



Human CDw293-BMPR1B recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6150
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Lys14-Arg126
TAG	C-Terminal Fc Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	41.34 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	O00238
ALIAS/SYNONYMS	BMP type-1B receptor, BMPR1B, CDw293, BMPR-1B, Bone morphogenetic protein receptor type-1B

