



Human CD172b-SIRPB1 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6151
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Ala369
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	41.08 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	O00241
ALIAS/SYNONYMS	SIRP-beta-1, Signal-regulatory protein beta-1, CD172 antigen-like family member B, SIRPB1, CD172b

