



Human CD277-BTN3A1 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6152
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Gly254
TAG	C-Terminal Fc Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	53.88 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	O00481
ALIAS/SYNONYMS	BTF5, BTN3A1, Butyrophilin subfamily 3 member A1, CD277

