



Human CD327-SIGLEC6 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6155
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Gln27-Val331
TAG	C-Terminal Fc Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	62.78 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	O43699
ALIAS/SYNONYMS	SIGLEC6, CD33L1, CD327, CD33L, Siglec-6, CDw327, CD33 antigen-like 1, OBBP1, Obesity-binding protein 1, Sialic acid-binding Ig-like lectin 6, OB-BP1

