



Human CD16b-FCGR3B recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6156
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Ser200
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	23.59 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	O75015
ALIAS/SYNONYMS	IGFR3, FCGR3B, CD16b, Low affinity immunoglobulin gamma Fc region receptor III-B, Fc-gamma RIII, IgG Fc receptor III-1, FCGR3, CD16B, Fc-gamma RIII-beta, Fc-gamma RIIIb, FCG3, FcRIII, FcR-10, FcRIIIb

