



## Human CD108-SEMA7A recombinant protein

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-P6159
<b>ORIGIN SPECIES</b>	Homo sapiens (Human)
<b>FRAGMENT</b>	His47-Ala648
<b>TAG</b>	C-Terminal Fc Tag
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Eukaryotic expression
<b>HOST</b>	Mammalian cells
<b>MOLECULAR WEIGHT</b>	96.6 kDa
<b>FORM</b>	
<b>BUFFER</b>	0.01M PBS, pH 7.4.
<b>STABILITY &amp; STORAGE</b>	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
<b>UNIPROT</b>	O75326
<b>ALIAS/SYNONYMS</b>	CDw108, Semaphorin-L, Sema L, SEMA7A, CD108, JMH blood group antigen, Semaphorin-K1, SEMAL, Sema K1, Semaphorin-7A, John-Milton-Hargen human blood group Ag

