



Human CD358-TNFRSF21 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6160
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Gln42-His349
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	36.66 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	O75509
ALIAS/SYNONYMS	DR6, Death receptor 6, Tumor necrosis factor receptor superfamily member 21, TNFRSF21, CD358

