



Human CD235a-GYPA recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6164
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Glu91
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	10.92 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P02724
ALIAS/SYNONYMS	PAS-2, Glycophorin-A, Sialoglycoprotein alpha, CD235a, MN sialoglycoprotein, GYPA, GPA

