



Human CD1a & B2M recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6166
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	CD1a[Met1-Val300] & B2M[Met1-Met119]
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	48.66 kDa(
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P06126 & P61769
ALIAS/SYNONYMS	T-cell surface glycoprotein CD1a, CD1A, T-cell surface antigen T6/Leu-6, CD1a, hTa1 thymocyte antigen, Beta-2-microglobulin, B2M

