



Human CD55 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6167
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Ser353
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	39.51 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P08174
ALIAS/SYNONYMS	Complement decay-accelerating factor, CR, DAF, CD55

