



Human CD66f-PSG1 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6170
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Pro419
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	48.19 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P11464
ALIAS/SYNONYMS	Pregnancy-specific glycoprotein 1, PSBG-1, Fetal liver non-specific cross-reactive antigen 1/2, B1G1, PSG1, PSGGA, PSG95, Pregnancy-specific beta-1 glycoprotein C/D, CD66 antigen-like family member F, PS-beta-C/D, CD66f, PS-beta-G-1, Pregnancy-specific be

