



Human CD222-IGF2R recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6171
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Asp1365
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	151.5 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P11717
ALIAS/SYNONYMS	IGF2R, M6P/IGF2 receptor, 300 kDa mannose 6-phosphate receptor, IGF-II receptor, M6PR, CI-MPR, CD222, Insulin-like Growth Factor proteins 2 receptor, MPR 300, M6P/IGF2R, MPRI, CI Man-6-P receptor, Insulin-like Growth Factor proteins II receptor, Cation-in

