



Human CD42c-GP1BB recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6172
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Cys147
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	16.37 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P13224
ALIAS/SYNONYMS	Antigen CD42b-beta, GP-Ib beta, GPIbB, GPIb-beta, GP1BB, CD42c, Platelet glycoprotein Ib beta chain

