



Human CD99 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6174
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Asp122
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	13.18 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P14209
ALIAS/SYNONYMS	MIC2Y, E2 antigen, MIC2X, 12E7, Protein MIC2, CD99 antigen, MIC2, T-cell surface glycoprotein E2, CD99

