



## Human CD75-ST6GAL1 recombinant protein

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-P6176
<b>ORIGIN SPECIES</b>	Homo sapiens (Human)
<b>FRAGMENT</b>	Glu44-Cys406
<b>TAG</b>	C-Terminal His Tag
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Eukaryotic expression
<b>HOST</b>	Mammalian cells
<b>MOLECULAR WEIGHT</b>	44.41 kDa
<b>FORM</b>	
<b>BUFFER</b>	0.01M PBS, pH 7.4.
<b>STABILITY &amp; STORAGE</b>	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
<b>UNIPROT</b>	P15907
<b>ALIAS/SYNONYMS</b>	Alpha 2,6-ST 1, B-cell antigen CD75, CMP-N-acetylneuraminate-beta-galactosamide-alpha-2,6-sialyltransferase 1, SIAT1, Sialyltransferase 1, ST6GAL1, ST6Gal I, Beta-galactoside alpha-2,6-sialyltransferase 1, ST6Gall

