



## Human ACVR2A recombinant protein

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-P6182
<b>ORIGIN SPECIES</b>	Homo sapiens (Human)
<b>FRAGMENT</b>	Met1-Pro134
<b>TAG</b>	C-Terminal Fc Tag
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Eukaryotic expression
<b>HOST</b>	Mammalian cells
<b>MOLECULAR WEIGHT</b>	41.54 kDa
<b>FORM</b>	
<b>BUFFER</b>	0.01M PBS, pH 7.4.
<b>STABILITY &amp; STORAGE</b>	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
<b>UNIPROT</b>	P27037
<b>ALIAS/SYNONYMS</b>	ACVR2, ACTRIIA, ACTR-IIA, Activin receptor type IIA, ACVR2A, Activin receptor type-2A

