



Human CD32b-FCGR2B recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6184
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Thr43-Pro217
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	22.51 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P31994
ALIAS/SYNONYMS	IGFR2, FCG2, Fc-gamma RII-b, CDw32, FcRII-b, IgG Fc receptor II-b, CD32, Fc-gamma-RIIb, Low affinity immunoglobulin gamma Fc region receptor II-b, FCGR2B

