



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Human CD118-LIFR recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6189
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Ser833
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	95.74 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P42702
ALIAS/SYNONYMS	LIFR, LIF-R, CD118, LIF receptor, Leukemia inhibitory factor receptor



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.