



Human CD151 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6192
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Ala113-Arg221
TAG	N-Terminal Fc Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	40.3 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P48509
ALIAS/SYNONYMS	CD151, PETA-3, Membrane glycoprotein SFA-1, GP27, Tspan-24, Platelet-endothelial tetraspan antigen 3, CD151 antigen, Tetraspanin-24, TSPAN24

