



Human CD239-BCAM recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6193
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Ala547
TAG	C-Terminal Fc Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	85.79 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P50895
ALIAS/SYNONYMS	Basal cell adhesion molecule, LU, MSK19, Auberger B antigen, CD239, Lutheran blood group glycoprotein, Lutheran antigen, F8/G253 antigen, B-CAM cell surface glycoprotein, BCAM

