



Human CD296-ART1 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6194
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Ser295
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	33.85 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P52961
ALIAS/SYNONYMS	CD296, ARTC1, GPI-linked NAD(P)(+)-arginine ADP-ribosyltransferase 1, ADP-ribosyltransferase C2 and C3 toxin-like 1, Mono(ADP-ribosyl)transferase 1, ART1

