



## Human CD85g-LILRA4 recombinant protein

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-P6197
<b>ORIGIN SPECIES</b>	Homo sapiens (Human)
<b>FRAGMENT</b>	Met1-Asn446
<b>TAG</b>	C-Terminal His Tag
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Eukaryotic expression
<b>HOST</b>	Mammalian cells
<b>MOLECULAR WEIGHT</b>	50.13 kDa
<b>FORM</b>	
<b>BUFFER</b>	0.01M PBS, pH 7.4.
<b>STABILITY &amp; STORAGE</b>	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
<b>UNIPROT</b>	P59901
<b>ALIAS/SYNONYMS</b>	CD85g, ILT7, Immunoglobulin-like transcript 7, LILRA4, Leukocyte immunoglobulin-like receptor subfamily A member 4, ILT-7, CD85 antigen-like family member G

