



Human CD156b-ADAM17 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6198
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Arg215-Asn671
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Insect
MOLECULAR WEIGHT	54.14 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P78536
ALIAS/SYNONYMS	TACE, Snake venom-like protease, ADAM17, CD156b, CSVP, ADAM 17, Disintegrin and metalloproteinase domain-containing protein 17, TNF-alpha convertase, TNF-alpha-converting enzyme

