



## Human CD167a-DDR1 recombinant protein

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-P6201
<b>ORIGIN SPECIES</b>	Homo sapiens (Human)
<b>FRAGMENT</b>	Met1-Thr416
<b>TAG</b>	C-Terminal His Tag
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Eukaryotic expression
<b>HOST</b>	Mammalian cells
<b>MOLECULAR WEIGHT</b>	46.87 kDa
<b>FORM</b>	
<b>BUFFER</b>	0.01M PBS, pH 7.4.
<b>STABILITY &amp; STORAGE</b>	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
<b>UNIPROT</b>	Q08345
<b>ALIAS/SYNONYMS</b>	Protein-tyrosine kinase RTK-6, Protein-tyrosine kinase 3A, NTRK4, EDDR1, CAK, Tyrosine kinase DDR, Mammary carcinoma kinase 10, TRKE, Epithelial discoidin domain-containing receptor 1, CD167a, TRK E, Epithelial discoidin domain receptor 1, HGK2, CD167 ant

