



Human CD215-IL15RA recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6203
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Ile31-Thr205
TAG	C-Terminal Fc Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	48.02 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q13261
ALIAS/SYNONYMS	IL-15 receptor subunit alpha, sIL-15R-alpha, IL-15RA, CD215, sIL-15 receptor subunit alpha, Interleukin-15 receptor subunit alpha, IL-15R-alpha, sIL-15RA, IL15RA

