



Human CD150-SLAMF1 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6204
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Lys236
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	27.31 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q13291
ALIAS/SYNONYMS	SLAM family member 1, IPO-3, CDw150, SLAMF1, SLAM, CD150, Signaling lymphocytic activation molecule

